# remel

## BactiCard® Candida

### INTENDED USE

Remel BactiCard® Candida is recommended for use in qualitative procedures to aid in the identification of *Candida albicans* by the rapid detection of L-proline aminopeptidase and  $\beta$ -galactosaminidase.

### SUMMARY AND EXPLANATION

*Candida albicans* is the species of yeast most frequently isolated from clinical specimens. Traditionally, *C. albicans* has been presumptively identified by the germ tube test.<sup>1</sup> Studies have indicated that up to 5% of clinical isolates of *C. albicans* are germ tube negative. Additionally, early pseudohyphae of *Candida tropicalis* may be confused with germ tube formation. Perry and Miller evaluated the use of 4-methylumbelliferyl-N-acetyl- $\beta$ -D-galactosaminide to detect an enzyme produced by *C. albicans*,  $\beta$ -galactosaminidase (MUGAL).<sup>2</sup> In further testing, Perry et al. modified the test substrate and added a second substrate for detection of L-proline aminopeptidase (PRO).<sup>3</sup> The combined use of two substrates provided for a method with improved sensitivity and specificity to aid in the identification of *Candida* spp. While other species of yeast were found to be positive for one enzyme or the other, *C. albicans* tested positive for both enzymes. This substrate combination provides a rapid, nonmicroscopic method to aid in the identification of *C. albicans* as effectively as the germ tube test.<sup>4</sup>

### PRINCIPLE

BactiCard® Candida consists of two separate test circles. The MUGAL circle contains the substrate, 4-methylumbelliferyl-N-acetyl- $\beta$ -D-galactosaminide. In the presence of the enzyme,  $\beta$ -galactosaminidase, this substrate releases 4-methylumbelliferone, a highly fluorescent compound detectable by longwave ultraviolet light. The PRO circle contains the substrate, L-proline- $\beta$ -naphthylamide. In the presence of the enzyme, L-proline aminopeptidase, the substrate is hydrolyzed and reacts with the Color Developer to form a red color.

### REAGENTS AND MATERIALS SUPPLIED

1. **BactiCard® Candida Test Cards** (25)  
Reactive Ingredients:  
L-Proline- $\beta$ -naphthylamide (PRO)  
4-Methylumbelliferyl-N-acetyl- $\beta$ -D-galactosaminide (MUGAL)
2. **Rehydrating Fluid** (1 bottle, 6.0 ml) with preservative
3. **Color Developer** (1 bottle, 6.0 ml)  
Reactive Ingredient: 0.02% p-Dimethylaminocinnamaldehyde
4. **MUGAL Reagent** (1 bottle, 6.0 ml)  
Reactive Ingredient: 4% Sodium Hydroxide
5. **Applicator Sticks** (1 vial)
6. **Instructions for Use** (IFU)

### PRECAUTIONS

This product is for *In Vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully.

### Caution!

1. MUGAL Reagent may cause irritation to skin, eyes, and respiratory system.
2. Color Developer is toxic and harmful by inhalation, contact with skin, or if swallowed. May impair fertility or cause harm to unborn child.
3. Refer to Material Safety Data Sheet for detailed information on reagent chemicals.

### STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 2-8°C until used. Do not freeze or overheat. Allow product to equilibrate to room temperature before use. Do not incubate prior to use. Once opened, maintain

unused test cards at 2-8°C in the foil pouch provided. Seal to protect from moisture.

### PRODUCT DETERIORATION

This product should not be used if (1) the color has changed, (2) the expiration date has passed, or (3) there are other signs of deterioration.

**Note:** Precipitate formation in the Color Developer at storage temperatures is common and does not affect performance of the test.

### SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.<sup>5</sup>

### MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Longwave ultraviolet light (Woods Lamp).

### PROCEDURE

Test isolates acceptable for testing with BactiCard® Candida are 18-72 hour, pure cultures of yeast.

1. Label BactiCard® Candida with the specimen identification.
2. Inoculate only one isolate per test card.
3. Add one (1) drop of BactiCard® Candida Rehydrating Fluid to each test circle. Do not oversaturate the test area.
4. Inoculate each test circle with a visible inoculum of the test isolate using an applicator stick (provided).
5. Incubate at room temperature for five (5) minutes.
6. Add one (1) drop of BactiCard® Candida Color Developer to the PRO test circle.
7. Observe for a red color development within 30 seconds, indicating a positive test.
8. Add one (1) drop of BactiCard® Candida MUGAL Reagent to the MUGAL test circle.
9. Observe the MUGAL test circle in a darkened room with a longwave ultraviolet light for bright blue fluorescence, indicating a positive test.

### INTERPRETATION

#### PRO Test:

Positive Test - Red color development within 30 seconds

Negative Test - No color change

#### MUGAL Test:

Positive Test - Bright blue fluorescence

Negative Test - No fluorescence

### EXPECTED VALUES

*C. albicans* is positive for both tests. If either the PRO or the MUGAL is negative or if both tests are negative, further testing is required for identification of the yeast isolate.

**Note:** *Candida dubliniensis* is also positive for both PRO and MUGAL. Growth at elevated temperatures and colony morphology on differential media have been reported to facilitate differentiation of *C. albicans* and *C. dubliniensis*.<sup>6,7</sup>

### QUALITY CONTROL

All lot numbers of BactiCard® Candida have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL	INCUBATION	RESULTS	
		PRO	MUGAL
<i>Candida albicans</i> ATCC® 10231	Aerobic, 5 min. @ Room Temperature	+	+
<i>Cryptococcus neoformans</i> ATCC® 32045	Aerobic, 5 min. @ Room Temperature	-	-

## ENGLISH

### PERFORMANCE CHARACTERISTICS

BactiCard® Candida was evaluated with 711 strains of yeast and demonstrated greater than 99% agreement with conventional identification methods. All 152 strains of *C. albicans* tested were PRO and MUGAL positive. Five of 152 strains were germ tube negative (3.3%). Of the 559 additional yeast isolates evaluated, six strains were positive for both enzymes (3 of 13 *Candida rugosa*, 2 of 13 *Candida stellatoidea*, and 1 of 6 *Trichosporon beigeli*, which is morphologically distinct from *C. albicans*).

### LIMITATIONS

1. This test is only part of the overall scheme for identification of yeasts. Further testing is required for definitive identification. Consult appropriate references for further instructions.<sup>1,5</sup>
2. Testing cultures mixed with bacteria or other species of yeast may yield aberrant results.
3. Testing cultures less than 18 hours old may result in false-negative reactions.
4. Extending the incubation of the test card beyond the specified time parameter may result in a false-positive reaction.
5. Some species of yeast may appear to fluoresce yellow on the MUGAL test circle and should not be confused with the blue fluorescence of a true positive reaction.
6. *C. dubliniensis*, first described in 1995, also produces germ tubes and chlamydospores.<sup>8</sup> Growth at elevated temperatures and morphology on differential media have been found to facilitate differentiation of *C. albicans* and *C. dubliniensis*.<sup>6,7</sup>

### BIBLIOGRAPHY

1. Winn, W., S. Allen, J. William, E. Koneman, G. Procop, P. Schreckenberger, and G. Woods. 2006. Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6<sup>th</sup> ed. Lippincott Williams and Wilkins, Baltimore, MD.
2. Perry, J.L. and G.R. Miller. 1987. J. Clin. Microbiol. 25:2424-2425.
3. Perry, J.L., G.R. Miller, and D.L. Carr. 1990. J. Clin. Microbiol. 28:614-615.
4. Dealler, S.F. 1991. J. Clin. Microbiol. 29:1081-1082.

5. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller. 2007. Manual of Clinical Microbiology. 9<sup>th</sup> ed. ASM Press, Washington, D.C.
6. Al Mosaid, A., D. Sullivan, I.F. Salkin, D. Shanley, and D.C. Coleman. 2001. J. Clin. Microbiol. 39:323-327.
7. Wabale, V.R., A.S. Kagal, R.S. Mani, and R. Bharadwaj. 2007. Indian J. Med. Microbiol. 25:304-305. Retrieved October 1, 2008 from: <http://www.ijmm.org/text.asp?2007/25/3/304/34787>.
8. Sullivan, D.J., T.J. Westerneng, K.A. Haynes, D.E. Bennett, and D.C. Coleman. 1995. Microbiology. 141:1507-1521.

### PACKAGING

REF R21106, BactiCard® Candida ..... 25 Tests/Kit

### Symbol Legend

REF	Catalog Number
IVD	In Vitro Diagnostic Medical Device
LAB	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
EC REP	European Authorized Representative

BactiCard® is a registered trademark of Remel Inc.  
ATCC® is a registered trademark of American Type Culture Collection.

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BactiCard® es una marca registrada de Remel Inc.